

Pipetting Techniques for Accurate Results



What is Air Displacement vs. Positive Displacement Pipetting?

An air displacement pipette, also known as an air cushion pipette, is a piston-operated instrument used to accurately aspirate and dispense precise liquid volumes. The piston and liquid in the plastic tip are separated by an air cushion allowing for contamination-free liquid handling. As the piston is moved up and down, a negative or positive pressure is created. This movement results in liquid being aspirated into or expelled from the tip.

Air displacement pipettes are ideal for liquids with similar density, viscosity, and vapor pressure (volatility) to water. For slightly volatile, viscous, or foaming liquids we recommend using the reverse pipetting technique. For highly volatile, viscous, or foaming liquids we recommend using a positive displacement instrument.

FORWARD PIPETTING



Aspirating the Sample

- Set the volume to be aspirated and on the pipette.
- Press the pipetting button to the first stop.
- Hold the pipette vertically and immerse the tip 1-4mm below the liquid surface. For larger volume pipettes (>1000µL) use an immersion depth of 3-6mm.
- Slowly release the pipetting button back to its original position.
- With tip immersed in liquid, wait several seconds until set volume is completely aspirated.

The ISO8655 standard requires pre-wetting the tip once prior to performing the pipetting procedure [additional technical information on ISO8655 calibration standards](#)

Dispensing the Sample

- Hold the pipette at a 30-45° angle touching the wall of the container.
- Push the pipetting button to the first stop and hold.
- Push the pipetting button to the second stop (blow-out) while wiping the tip against the wall.
- Remove the pipette tip from the container, and then release the pipetting button back to its original position.



PIPETTING TIPS

Attaching Pipette Tips

Proper tip fit will ensure that only a modest force is required to attach a tip to the pipette. "Tapping" and "rocking" the pipette to attach tips suggests an improper tip fit and can result in the pipette needing to be adjusted and recalibrated.

Environment Matters

Environmental conditions (e.g. humidity, temperature, etc.) can impact the properties of the air cushion leading to inaccuracies. We recommend keeping environmental conditions consistent while pipetting.

Use Quality Pipette Tips

High quality tips provide an airtight seal without the need for excessive force (tapping or rocking) to attach the tip to the pipette. Consider tips made of superior materials, that are free of additives, and manufactured from molds that avoid departing agents or the use of molding aids.

Ejection Force

Thin-walled, flexible pipette tips are designed to fit the pipette securely while allowing for minimal attachment and ejection forces. Low ejection force improves ergonomics which decreases stress in the user.


Preparing the Pipette

Pipetting Accuracy is best at nominal (maximum) volume, so for optimal results, select a pipette with nominal volume closest to aspirated liquid volume. Ensure the volume setting is correct and attach a BRAND® or other high quality pipette tip.

REVERSE PIPETTING

Aspirating the Sample

- Pre-wet the tip 3-5 times by slowly aspirating and dispensing liquid in the tip (forward pipetting).
- Press the pipetting button to the second stop.
- Hold the pipette vertically and immerse the tip 1-4mm below the liquid surface. For larger volume pipettes (>1000uL) use an immersion depth of 3-6mm.
- Slowly release the pipetting button back to its original position, paying attention to avoid bubbles in the tip.
- With tip immersed in liquid, wait several seconds until set volume is completely aspirated.



Dispensing the Sample

- Hold the pipette at a 30-45° angle touching the wall of the container.
- Push the pipetting button to the first stop and hold for several seconds.
- Remove the pipette tip from the container while holding pipetting button at the first stop. (Note: Excess liquid will remain in tip and should not be dispensed)
- Excess liquid can be returned to original solution by pressing pipetting button to second stop or ejected with the tip.




BRAND® QUALITY


[BRAND Transferpette® S](#) pipettes and [BRAND tips](#) are the ideal combination for quality and accuracy in pipetting. The precision workmanship, simple operation, and ergonomic comfort of the Transferpette S make it the performance standard in pipetting.




ACCURACY AND PRECISION




Inaccurate and imprecise



Accurate but imprecise



Inaccurate but precise



Accurate and precise

- Calibration is the process of determining the performance of an instrument without mechanical or electronic adjustment.
- Adjustment is the manipulation of the instrument post-calibration so that the instrument is aligned within the specified tolerances.
- Accuracy of measurement is how close the result comes to the target value.
- Precision (reproducibility) describes how closely grouped results are in a set of measurements, in units of volume.