

Reliable Tools for Vape & E-liquid Labs

The Standard in Bottletop Dispensers

Dispensette® S dispensers

- Rugged: no seals to wear out!
- Compatible with most reagents, including nicotine
- Simple to clean & calibrate
- Recirculation valve reduces splashing and helps conserve reagents
- Autoclavable to 250 °F (121 °C).
- Adapters for common bottle threads included
- Optional flexible discharge tube simplifies serial dispensing



Dispensette S	No recirculation valve	Recirculation valve
	Volume, mL	Cat. No.
Dispensette S Digital		
0.1-1	4600310	4600311
0.2-2	4600320	4600321
0.5-5	4600330	4600331
1-10	4600340	4600341
2.5-25	4600350	4600351
5-50	4600360	4600361
Dispensette S, Analog-adjustable		
0.1-1	4600100	4600101
0.2-2	4600120	4600121
0.5-5	4600130	4600131
1-10	4600140	4600141
2.5-25	4600150	4600151
5-50	4600160	4600161
10-100	4600170	4600171

Fixed volume models available. For strong acids and nonpolar organics, ask for the BRAND Dispensette® S Organic dispensers.

Flexible discharge tube

	Cat. No.
For Nominal Volume, mL. PTFE, coiled, length 800 mm (= 31.5") with safety handle.	
For Dispensette® S and Dispensette® S Organic	
2, 5, 10	708132
25, 50, 100	708134
For seripettor	
2, 10	704522

*Not suitable for use with peroxides due to platinum-iridium valve spring.

A Rugged, Economical Tool for Routine Dispensing

seripettor®

- Dispenses most non-aggressive liquids such as propylene glycol and vegetable glycerin
- Convenient cleaning and maintenance
- Notched volume selector allows for quick and exact volume adjustments
- Comes with thread adapters to fit most common lab bottles



Volume, mL	Subdivision, mL	A* <±%	CV* ≤%	Cat. No.
seripettor				
0.2-2	0.04	1.2	0.2	4720120
1-10	0.2	1.2	0.2	4720140
2.5-25	0.5	1.2	0.2	4720150

* The value of accuracy and coefficient of variation are final test values referring to the delivered nominal volume, instrument and distilled water at equilibrium with ambient temperature (20 °C/68° F) and with smooth operation. A*=Accuracy, CV*=Coefficient of Variation

Serial dispensing

The optional flexible discharge tube with safety handle speeds serial dispensing tasks, and permits fast and precise dispensing even into narrow test tubes. Functions of recirculation valve and safety discharge system are fully maintained with the flexible discharge tube.






Ethylene-acrylate Coated Glass Amber Bottles

The ethylene-acrylate coating envelopes the glass bottle like a protective skin. Should the bottle break, the splintering hazards are reduced considerably.

Description	Cat. No.
Threaded Bottle, PP screw cap, LDPE pouring ring, 250 mL, 33 mm	704004
Threaded Bottle, PP screw cap, LDPE pouring ring, 500 mL, 33 mm	704006
Threaded Bottle, PP screw cap, LDPE pouring ring, 1000 mL, 45 mm	704008
Bottle Stand, PP, 225 mm x 160 mm x 325 mm	704275



HandyStep® Repeating Pipettes

HandyStep S

- Versatile, reliable mechanical repeating pipette
- Comfortable, lightweight and well-balanced for one-handed operation
- Flexible – accepts third party tips

HandyStep touch

- Electronic repeating pipette with simple touchscreen operation
- Easy - automatic tip size recognition of BRAND® PD-Tip //s
- Efficient – accepts third party tips

Description	Cat. No.
HandyStep S	705110
HandyStep touch	705200
HandyStep touch S	705210



BRAND® PD-tip™ // Syringe Tips

BRAND encoded PD-Tip // syringe tips are manufactured from high-quality virgin plastics for use with all standard repeating pipettes in a wide variety of applications.

- Compatible with most pipettes
- Reliable
- Manufactured from high quality materials
- Batch certificate of performance

Description	Cat. No.	
HandyStep S Repeating Pipette	705110	
HandyStep touch Repeating Pipette	705200	
PD-Tip // size		
0.1 mL (100 tips)	705700	Non-sterile tips
0.5 mL (100 tips)	705702	BIO-CERT® Pure*
1.0 mL (100 tips)	705704	705730
1.25 mL (100 tips)	705706	705732
2.5 mL (100 tips)	705708	705734
5 mL (100 tips)	705710	705736
10.0 mL (100 tips)	705712	705738
12.5 mL (100 tips)	705714	705740
25 mL (50 tips)	705716	705742
25 mL (25 tips)	—	705744
50 mL (25 tips)	705718	—
Multi-pack (20 tips each 0.5, 1, 1.25, 2.5, 5, 10 & 12.5 mL)	705720	705746
		705748

*BIO-CERT is certified sterile, and DNA-, RNase-, ATP- and endotoxin-free.